(4% icodextrin solution)

MANUFACTURER:

Baxter AG
A-1221 Vienna, Austria

DESCRIPTION

ADEPT® is a single use, sterile, clear, colourless to pale yellow fluid for intraperitoneal administration containing icodextrin at a concentration of 4% w/v in an electrolyte solution.

Each 1 litre of solution contains:

Icodextrin	40 g
Sodium Chloride	5.4 g
Sodium Lactate	4.5 g
Calcium Chloride	257 mg
Magnesium Chloride	51 mg
Theoretical osmolarity 278 milliosmoles per litre	

Ionic composition (approximately) per litre:

Sodium	133 mmol
Calcium	1.75 mmol
Magnesium	0.25 mmol
Chloride	96 mmol
Lactate	40 mmol

ADEPT is packaged in flexible polyvinylchloride bags containing either 1 litre or 1.5 litres of solution.

INDICATIONS

ADEPT is intended for use as an intraperitoneal instillate for the reduction of adhesions following abdominal surgery, and should be used as the irrigant during the course of that surgery.

ACTIONS

ADEPT performs its function through a physical effect by providing a temporary separation of peritoneal surfaces by hydroflotation. This minimises tissue apposition during the critical period of fibrin formation and mesothelial regeneration following surgery, thereby providing a barrier to adhesion formation.

Icodextrin is an α -1,4-linked glucose polymer which, when administered intraperitoneally as a 4% solution, is capable of maintaining a reservoir of fluid within the peritoneal cavity for up to 3–4 days. A 7.5% solution of icodextrin has been used extensively on a daily basis as a peritoneal dialysis solution for the treatment of chronic renal failure.

When given intraperitoneally, the polymer is largely retained within the peritoneal cavity. Some absorption occurs from the peritoneum into the systemic circulation where it is metabolised by amylase to smaller oligosaccharides, ultimately maltose and by maltase to glucose.

ADEPT has been shown to reduce significantly the incidence, extent and severity of post surgical adhesions in animal models (rabbit double uterine horn and rabbit sidewall models) when used as a lavage during surgery and as an instillate post-operatively.

CONTRAINDICATIONS

ADEPT should not be used in patients with a known allergy to starch based polymers or in patients with maltose or isomaltose intolerance.

PRECAUTIONS**Use in Children**

ADEPT is not recommended for use in children.

Pregnancy and Lactation

There are limited data available from animal studies on the effects of icodextrin on reproduction or lactation and therefore ADEPT should not be used during pregnancy or lactation.

Women of childbearing potential should be treated with ADEPT only when adequate contraceptive precautions have been taken.

Interactions with Other Medicaments

The primary intended function of ADEPT is not to administer medicinal products. However, the bag has an injection port, which may be used for administration of drugs, if required.

A range of antibiotics, including vancomycin, cephazolin, ampicillin, flucloxacillin, ceftazidime, gentamycin and amphotericin, have shown no evidence of incompatibility with ADEPT.

UNDESIRABLE EFFECTS

Undesirable effects are those typically seen following surgery. In patients receiving icodextrin 7.5% solution as part of a peritoneal dialysis regimen and on multi therapy, there have been common reports of skin reactions, including rash and pruritus. Occasionally these rashes have been associated with exfoliation.

There have been rare reports of hypersensitivity reactions in patients treated with ADEPT.

There have been rare reports of vulval oedema following the administration of ADEPT. The reaction generally resolves spontaneously within a few days. Oedema is a recognised event associated with the use of fluids for irrigation and instillation in laparoscopic surgery.

DIRECTIONS FOR USE

ADEPT is administered into the peritoneal cavity during abdominal surgery, being used as an irrigant solution during the course of surgery. Once the surgeon has completed the surgical procedure(s) and removed all packs and sponges, the cavity is aspirated of all remaining fluid. A final volume of at least 1 litre of ADEPT is then introduced into the cavity before closure of the cavity/removal of the scope.

ADEPT should be warmed to approximately body temperature prior to use, using a device specifically intended for warming solutions in operating theatres. ADEPT can be kept in a warmer at 37°C for up to 14 days, provided it is not removed and then replaced. At all other times, storage below 4°C or above 30°C is not recommended.

Using standard operating room technique:

1. Remove the outer wrap from the ADEPT bag and hang the sterile bag of solution on a stand.
2. Remove the twist-off tab from the spike port and insert a standard giving set for connection to a laparoscope or a giving set for dispensing the solution directly into the abdominal cavity in the case of laparotomy.
3. ADEPT should be used intra-operatively as an irrigant solution, and as a post-operative instillate. The solution will flow through a giving set (and through laparoscopes), or it can be dispensed into a sterile basin and applied using a syringe and cannula.
4. When used as an intra-operative irrigant solution, at least 100 mls of ADEPT should be introduced to the cavity every 30 minutes.
5. Remove remaining fluid before introducing the final instillation.
6. For the final instillation of ADEPT, prior to closure of the abdominal cavity or removal of the laparoscope, at least one litre (a new bag of ADEPT if 1 litre bags are being used) should be used. Direct the solution at the operative sites in the first instance, the remainder being distributed throughout the cavity.
7. Dispose of the bag and any unused portion of the solution following normal operating room biological hazard procedures.

PRECAUTIONS FOR USE

ADEPT must be used as directed by a physician. It must not be used unless the solution is clear and the container undamaged.

Any unused portion of solution should be discarded. ADEPT is not to be used for intravenous infusion.

STORAGE

ADEPT should not be stored above 30°C. Do not refrigerate or freeze.

PRESENTATION

ADEPT is packaged in single use, flexible polyvinylchloride bags, fitted with connecting ports, containing 1 litre or 1.5 litres of solution. The product is presented sterile (by heating in an autoclave). The bags are packaged in cartons of 10 x 1 litre or 5 x 1.5 litres.

ADEPT is a Registered Trade Mark of Innovata plc

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