VERSION 5

0471/02/1

Baxter

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Vertrieb
Baxter AG
CH-8604 Volketswil

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STERILE EO Sterilisationsmethode: Ethylenoxid

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INSTRUCTIONS FOR USE

GentaFleece collagen fleece
Collagen sponge with antibiotic protection for surgical use – haemostatic, sterile –

Composition
1 sponge 6 x 4 x 0,5 cm contains:
Collagen from equine tendons 67.2 mg
Gentamicin sulphate 48 mg corresponding to 26.4 – 34.32 mg gentamicin

1 sponge 12 x 8 x 0,5 cm contains:
Collagen from equine tendons 268.8 mg
Gentamicin sulphate 192 mg corresponding to 105.6 – 137.28 mg gentamicin

1 cm² sponge 0.5 cm in thickness contains:
Collagen from equine tendons 2.8 mg
Gentamicin sulphate 2 mg corresponding to 1.10 – 1.43 mg gentamicin

Pharmaceutical form
Absorbable implant

Indications

GentaFleece collagen fleece is indicated for local haemostasis of capillary bleeding from parenchymatous tissue in regions with a high risk of infection. GentaFleece collagen fleece may be applied to defect cavities and other residual cavities during surgery, as for example in the sacral cavity after rectum amputation, or for soft-tissue abscesses resulting from pilonidal sinuses.

GentaFleece collagen fleece may be applied for haemostasis in clean, clean-contaminated and contaminated wound cavities in the presence of diffuse, capillary, arterial-venous, or arterial or venous bleeding, extensive capillary bleeding from parenchymatous organs, or as a supportive measure for other procedures for haemostasis.

GentaFleece collagen fleece may also be applied in combination with fibrin adhesives.

Dosage and administration
After removing the focus of inflammation GentaFleece collagen fleece is applied into the defect cavity as required. GentaFleece collagen fleece is loosely inserted as a flat sponge – either rolled or folded – or, in cases of bone defects, applied mixed with bone chips.

The collagen sponge is absorbable and may not be removed. Haemostasis occurs on a physical basis. The contact of collagen with blood leads to aggregation of platelets which precipitate in large numbers on the collagen matrix, disintegrate, and release clotting factors which, together with plasma factors, facilitate the formation of fibrin. The collagen matrix gives additional reinforcement to the blood clot.

Due to its structure GentaFleece collagen fleece is capable of absorbing large quantities of fluid. During this purely mechanical process of absorption of secretion, rejected materials such as bacteria and fibrin clots are also removed. This accelerates the formation of granulation tissue. The size of the defect and the weight of the patient determine the dosage of GentaFleece collagen fleece. Subject to these parameters, 1 to 3 sponges (10 x 10 cm) are generally applied for patients up to 50 kg, or a maximum of 5 sponges (10 x 10 cm) for patients over 50 kg body weight. The number of sponges is tailored according to an appropriate number of GentaFleece collagen fleece sponges of the size 5 x 5 cm may be used.

GentaFleece collagen fleece is applied in a dry state with light pressure to achieve better adhesion. Before use of the collagen fleece, the area to be treated should be disinfected. The release of the incorporated gentamicin occurs simultaneously with the dissolving of the collagen sponge. This protects the implant from external contamination by bacteria ascending along the drainage and from contamination by bacteria spread or introduced during the surgical procedure. In cases of diffuse bleeding (e.g. in the sacral cavity) it may prove useful to apply the implant by pressing it into the affected areas of the wound. Overpacking the wound with collagen may result in the formation of a seroma. Care should be taken to apply the sponges into the tissue of the wound as far as possible.

When using GentaFleece collagen fleece in combination with a fibrin adhesive, apply the adhesive to the surface of the sponge which will lie against the area being treated.

Contraindications
GentaFleece collagen fleece must not be used in patients with a known hypersensitivity to collagen or gentamicin or other aminoglycosides. GentaFleece collagen fleece should only be used as an adjunct to other procedures for haemostasis in the presence of bleeding requiring control by ligatures or in instances of larger arterial and/or venous haemorrhage requiring suture ligatures.

Warnings

GentaFleece collagen fleece should only be used with strict medical indications in patients with impaired renal function as well as in patients with autoimmune diseases or neuromuscular diseases such as Parkinson's or myasthenia gravis.

Although toxic serum levels are not reached during the use of GentaFleece collagen fleece, both plasma gentamicin levels and renal function should be monitored.

The concomitant systemic administration of aminoglycoside antibiotics should be avoided or these should only be given under strict monitoring of serum gentamicin levels and renal function. Cross-sensitivity to aminoglycoside antibiotics do exist.

There is no reported experience with the use of GentaFleece collagen fleece in patients with immunological or connective tissue diseases. Although it has not been proven that the use of animal collagen results in exacerbation of the disease, GentaFleece collagen fleece should only be used with strict medical indications in these patients.

There is no reported experience with the use of GentaFleece collagen fleece in children.

Special precautions for use
For the product, conventional procedures for the application of a haemostatic agent should be observed. For the risk-free use of GentaFleece collagen fleece the following precautions for use should be observed:

- Soaking or moistening GentaFleece collagen fleece prior to implantation may result both in a loss of efficacy of its haemostatic properties as well as in a loss of its self-protective property through premature precipitation of the mildly water-soluble gentamicin sulphate. A reduction in the sponge's blood absorption capacity following the soaking of GentaFleece collagen fleece results in a reduction of its physical properties of platelet aggregation and thus to a decisive interference with haemostasis.

- The desired platelet aggregation is also achieved by lining the wound with GentaFleece collagen fleece. The simultaneous insertion of a suction and irrigation drain is strictly prohibited as this would lead to contamination by bacteria spread or introduced during the surgery by the collagen sponge against infection. The drain openings may also become blocked by fragments of collagen.

- GentaFleece collagen fleece should not be used alone. An additional antibiotic should be administered systemically.

Interactions with other medicinal products

Although only very low plasma levels are reached after implantation of GentaFleece collagen fleece, the interaction of gentamicin with other medicinal products should be considered.

The concomitant administration of aminoglycosides with loop diuretics, e.g. furosemide or etacrynic acid, should be avoided as these diuretics themselves have an ototoxic effect. If diuretics are simultaneously administered intravenously, they may increase the concentration of gentamicin in serum or tissue and thus increase toxicity.

The simultaneous or consecutive or topical application of potentially neurotoxic and/or nephrotoxic substances (e.g. colistin, other aminoglycosides, streptomycin, cefaloridine, viomycin, polymyxin B or polymyxin E, may increase the toxicity of gentamicin.

The simultaneous topical use of beta-lactam antibiotics may result in altered activity (inactivation).

The neuromuscular blocking tendency may be increased by the simultaneous application of muscle relaxants, e.g. d-tubocurarine, suxamethonium or pancuronium, as well as by ether. A neuromuscular blockade can best be reversed by the administration of calcium salts. Should a neuromuscular blockade occur under the simultaneous application of suxamethonium in the presence of acquired or genetically determined cholinesterase deficiency, then artificial respiration will be necessary and cholinesterase must be administered.

Pregnancy and lactation
No data on the safe use of GentaFleece collagen fleece during pregnancy and lactation are available.

Antibiotic substances which cross the placenta and reach the fetus. There is therefore a risk of intrauterine fetal damage.

Since aminoglycoside antibiotics also pass into breast milk, GentaFleece collagen fleece should also not be used during lactation, or breastfeeding should be ceased.

Effects on ability to drive and operate machines
None known yet.

Adverse reactions
There may be increased secretion once absorption of GentaFleece collagen fleece commences in the tissue. Overflow drains should therefore be used whenever possible in cases of extensive infected cavities.

Nephrotoxicity and ototoxicity may be potentiated by the simultaneous systemic administration of aminoglycoside antibiotics in or patients with impaired renal function.

Overdose
There is no appreciable increase in serum gentamicin levels when the recommended amounts of GentaFleece collagen fleece are applied. There is no information available on doses above these levels.

Plasma gentamicin levels and renal function should always be monitored when doses other than those recommended are given.

Peritoneal dialysis or haemodialysis should be considered in cases of severe intoxication.

Incompatibilities
None known yet.

Storage
Do not store above 25°C. Protect from moisture.

Keep out of the reach of children.

Shelf life, sterility
The expiry date and the batch number are imprinted on the retail pack and on the cover of the blister pack. GentaFleece collagen fleece must not be used after the expiry date. GentaFleece collagen fleece sponges removed from the protective package must not be resterilised and should therefore be discarded.

The contents of unused, yet opened or damaged, packets must not be resterilised and should therefore be discarded. GentaFleece collagen fleece is intended for single use only.

Pack sizes
Pack with 5 sponges à 24 cm² (6 x 4 cm)

REF B224400099999

Pack with 1 sponge à 96 cm² (12 x 8 cm)

REF B224420099999

Pack with 5 sponges à 96 cm² (12 x 8 cm)

REF B224440099999

Manufacturer
RESORBA Wundversorgung GmbH + Co. KG

Am Flachmoor 16
D-90475 Nürnberg
C E 0125

Distributor
Baxter AG
Industriestrasse 67
A-1220 Wien

Distributor
Baxter AG
CH-8604 Volketswil

Symbols and explanations

STERILE EO Sterilisation method: ethylene oxide

⊗ Do not resterilise, do not reuse

⚠ Attention, see instructions for use

LOT Batch number

📅 Expiry date year/month

REF Product number

25°C Upper temperature limit

Date of preparation: 03/2002

Baxter is a trademark of Baxter International Inc.

GentaFleece is a trademark of Baxter International Inc. and Baxter Healthcare S.A.

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