# TISSUDURA

Collagen biomatrix for dura regeneration

## Instructions for Use, read carefully!

The information given in these Instructions for Use is regularly updated according to the latest findings. Please read this information carefully prior to each use of product.

## COMPOSITION

1 cm<sup>2</sup> TISSUDURA collagen biomatrix contains 5,6 mg of native collagen fibrils of equine origin.

## INDICATIONS

TISSUDURA is used for dura reconstruction as a temporary dura replacement and as a biomatrix for dura regeneration. TISSUDURA can be used together with fibrin sealant.

## MODE OF ACTION

Equine collagen fibrils in TISSUDURA serve as a matrix to support tissue repair. Therefore, TISSUDURA acts as a biomatrix for dura replacement and dura regeneration by triggering the formation of a new extracellular matrix and rapid migration of fibroblasts.

The non-porous structure of TISSUDURA is already primarily water-tight.

The contact of collagen with blood results in the aggregation of platelets. When platelets adhere to the collagen matrix they degrade and release coagulation factors, which in combination with other plasma factors, lead to the formation of fibrin. Therefore the use of tissue adhesives for sealing is not mandatory when TISSUDURA is applied.

## DOSAGE AND ADMINISTRATION

TISSUDURA can be cut to the desired size and shape.

It is recommended that the TISSUDURA implant overlaps the dura defect by approximately 1 cm.

Prior to application TISSUDURA can be soaked in sterile physiological saline solution (approx. 5 min.), which can contribute to improved handling properties.

TISSUDURA can be soaked with therapeutic solutions such as antibiotics (see also "Drug Interactions") prior to application.

In special situations TISSUDURA may be applied without prior soaking.

TISSUDURA should be applied in compliance with neurosurgical principles. If necessary, fibrin sealant can be used for fixation of TISSUDURA.

The TISSUDURA implant should be prepared and placed in such a way that its structural integrity is not disrupted. It should be fitted and fixed in place in such a way as to prevent displacement.

Normally, suturing is not required, but tensionless, atraumatic sutures may be used if desired. If suturing is used, sealing of the sutures with fibrin sealant is recommended.

## CONTRAINDICATIONS

The product must not be used in patients with a known hypersensitivity to products of equine origin.

The application of TISSUDURA onto infected areas is contraindicated.

TISSUDURA should not be used in applications where the implant is subjected to an increased mechanical load or a high degree of stretching.

TISSUDURA should not be used in procedures involving large defects, where adequate fixation is unlikely.

TISSUDURA should not be used in the presence of wound healing disorders.

There is no data available on the use during pregnancy and lactation. Therefore, any use of the product during pregnancy and lactation has to be preceded by a cautious risk-benefit assessment.

## UNDESIRABLE EFFECTS

The literature reports an occasional aggravation of a pre-existent infection when collagen is implanted. Due to the xenogenic origin of the collagen, hypersensitivity reactions may occur despite the known low antigenic property.

## DRUG INTERACTIONS

Halogen-releasing antiseptic agents as well as tannins and caustics should not be used in combination with TISSUDURA.

Due to the specific effect of collagen on the platelet aggregation and plasma coagulation, the hemostatic efficacy of TISSUDURA may be reduced by substances inhibiting aggregation and coagulation or in the presence of coagulation disorders.

## SHELF LIFE AND STORAGE

Do not store above 25°C. Keep dry.

Expiry date: See imprint on outer package. Do not use after the expiry date.

After opening the sterile single pack, do not resterilise the remaining content, but discard the unused product.

Keep out of the reach and sight of children.

## DIMENSIONS AND PACK SIZES

TISSUDURA, packed in sterile single packs, is available in the following dimensions:

Dimensions	Pack sizes	REF Number
2,5 x 2,5 cm	1 x	B2246000999999
2,5 x 10 cm	1 x	B2246200999999
5 x 5 cm	1 x	B2246400999999
5 x 10 cm	1 x	B2246600999999
10 x 10 cm	1 x	B2246800999999

## SYMBOLS AND EXPLANATIONS

	Sterile, sterilisation method: ethylene oxide
	Do not resterilise, do not reuse
	Attention, see instructions for use
	Batch number
	Expiry date year/month
	Product number
	Maximum storage temperature
	Keep dry
	Manufacturer

Date of preparation: 08/2004

BAXTER is a trademark of Baxter International Inc.  
TISSUDURA is a trademark of Baxter Healthcare S.A.

	<b>Distributor:</b>
Resorba Wundversorgung GmbH & Co. KG Am Flachmoor 16 D-90475 Nuremberg	Baxter Healthcare S.A. Hertistrasse 2 CH-8304 Wallisellen

**CE 1275** = CE-Mark and identification number of Notified Body.  
Product conforms to the essential requirements of the medical device directive 93/42 EEC.

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BAXTER is een handelsmerk van Baxter International Inc. TISSUDURA is een handelsmerk van Baxter Healthcare S.A.