Baxter

S

TissuFleece E

Kollagenplatta

Bipacksedel, läs noggrant!

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Collagen fleece

Package leaflet, read carefully !

The information given in this package leaflet is regularly updated. Please read instructions carefully prior to use.

ACTIVE SUBSTANCE

Collagen fibrils

ACTIVE INGREDIENTS

1 cm² contains 2.8 mg of native collagen fibrils of equine origin.

OTHER INGREDIENTS

None

MODE OF ACTION

The contact of collagen with blood results in the aggregation of thrombocytes which then adhere in great numbers to the collagen matrix, disintegrate, and release coagulation factors which, in combination with plasma factors, lead to the formation of fibrin.

The collagen matrix additionally enhances the coagulum.

Owing to its structure, TissuFleece E collagen fleece is capable of absorbing large quantities of liquid. By way of this merely mechanical process, cast off materials such as bacteria and fibrin layers are absorbed. Consequently, the formation of granulation tissue is accelerated.

THERAPEUTIC INDICATIONS

- In wound management and as a haemostyptic agent in venous, capillary, diffuse, and extensive bleedings from parenchymatous organs, such as lungs, liver, spleen, kidney, prostate.
- In combination with surgical tissue sealants, for the sealing and stabilisation of sutures with an inherent risk of insufficiency, e.g., in heart and vascular surgery.
- For haemostasis following tooth extraction.

- As temporary skin replacement for painless covering, e.g., of leg ulcers.
- For the covering of skin graft donor sites.
- For temporary covering of grade 2 and 3 burns.
- For the restoration of defects and haemostasis after resection of spongiosa or other orthopaedic interventions.

CONTRAINDICATION

The product must not be used in patients with hypersensitivity to products of equine origin. The application onto infected areas is an absolute contraindication. There is no data available on the use during pregnancy and lactation. Therefore, any use of the product during pregnancy and lactation has to be preceded by a cautious assessment of the indication.

UNDESIRABLE EFFECTS

The literature reports occasional events suggesting that potentially present infections may be aggravated. Due to the non human origin of the collagen, hypersensitivity reactions may occur despite the known low antigenic property.

DRUG INTERACTIONS

Chlorine-releasing antiseptic agents as well as tannins and caustics should not be used in combination with collagen.

Due to the specific effect of TissuFleece E collagen fleece on thrombocyte aggregation and plasma coagulation, attention has to be paid to the fact that its efficacy will be reduced through substances inhibiting aggregation and coagulation.

POSOLOGY AND MODE OF APPLICATION

Unless otherwise prescribed, TissuFleece E collagen fleece is cut to the desired size and placed onto the wound without prior treatment. When used topically, it is fixed by applying a bandage, when used internally it is fixed, in order to enhance the haemostatic effect, either by applying surgical tissue sealant, such as GRF sealant, fibrinogen-thrombin-F-XIII sealant, or, if required, by placing conventional sutures.

In the treatment of leg ulcers, TissuFleece E collagen fleece is renewed 2–3 times per week. TissuFleece E collagen fleece can be applied dry or after soaking in physiological saline solution or, e.g., antiseptic agents (see also Drug Interactions). In the event of special indications, TissuFleece E collagen fleece can be soaked with therapeutic solutions such as antibiotics prior to application.

Store out of the reach of children!

ATTENTION

Do not sterilise unused sheets removed from the immediate packaging.

Due to the affinity of collagen to bleeding surfaces, dry instruments and gloves have to be used when applying TissuFleece E collagen fleece.

SHELF LIFE AND STORAGE

Store in a dry place below 25°C.

Do not sterilise opened individual packages.

Expiry date: See imprint on outer package.

Do not use after the expiry date!

DIMENSIONS UND PACK SIZES

TissuFleece E collagen fleece, sterile, is available in the following dimensions and pack sizes:

Dimensions	Pack sizes	REF Number
3,6 x 1,8 cm	20 pieces	B2240700999999
7 x 3 cm	5 pieces	B2240500999999
9 x 7 cm	5 pieces	B2240300999999
12 x 9 cm	5 pieces	B2240100999999

SYMBOLS AND EXPLANATIONS

	Sterile, sterilization method: ethylene oxide
	Do not resterilise, do not reuse
	Attention, see instructions for use
	Batch number
	Expiry date year/month
	Product number
	Maximum storage temperature
	Manufacturer

Date of preparation: 02/2003

Baxter is a trademark of Baxter International Inc.

TissuFleece is a trademark of Baxter International Inc. and Baxter Healthcare S.A.

Manufacturer:	Resorba Wundversorgung GmbH & Co. KG Am Flachmoor 16 D-90475 Nuremberg	Distributor:	Baxter AG Industriestrasse 67 A-1220 Vienna
----------------------	--	---------------------	---

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Baxter é uma marca registada da Baxter International Inc.
TissuFleece é uma marca registada da Baxter International Inc. E Baxter Healthcare S.A.