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### INSTRUCTIONS FOR USE

### DESCRIPTION

TRICOS is a resorbable bioactive bone substitute, which is able to be replaced by newly formed bone over time. TRICOS is a biphasic calcium phosphate ceramic constituted of 60 % Hydroxyapatite (HA) and 40 % beta Tricalcium phosphate (β-TCP). 3,5 cm<sup>3</sup> of TRICOS granules are contained in an application device.

HA is chemically comparable to biological apatite crystals. Chemically pure HA is considered as the least bioactive material among the calcium phosphate ceramics *in vitro* and *in vivo*. β-TCP is recognized to have a higher bioactivity and dissolution rate than HA *in vitro* and *in vivo*.

It has been shown *in vitro* and *in vivo* that a mixture of 60 % HA and 40 % β-TCP represents the ideal composition for a biphasic calcium phosphate ceramic for use as a resorbable bone substitute. TRICOS is presented in the form of porous granules with a particle size of 1-2 mm. The porosity is required for the biological *in vivo* exchanges, particularly for bone ingrowth and mineralization. There are two types of porosity contained in TRICOS:

- The mechanical properties of the material should be taken into account when subjecting to stress and, if necessary, to stabilize (with a fixator or by osteosynthesis).
- The implantation in a revision surgical site with non-resorbable fragments of material (such as polyethylene ligament waste, carbon fibers...) is not recommended.
- Pre and post-operative Procedures are performed according to hospital routine.
- To preserve the porous structure, the precautions for use must be strictly followed.
- When mixed with TISSEEL/TISSUCOL Kit, the 4 IU/ml Thrombin concentration must be used.

#### SIDE-EFFECTS

No side effects associated with this type of biomaterials have been reported.

#### STERILIZATION / STORAGE

TRICOS is a class III medical device. TRICOS is supplied in sterile package. It has undergone beta radiation, which makes it suitable for surgical use. TRICOS should be stored at room temperature.

- ⊗ Do not re-use

<b>STERILE</b>	<b>R</b>
----------------	----------

<b>REF</b>	<b>B303500099999</b>
------------	----------------------

#### 0123

This product is to be handled or implanted by trained qualified persons having read these instructions for use.

<b>Manufacturer:</b> Biomatlante SARL	<b>Distributore:</b> Baxter Healthcare S.A.
Z.A. Les Quatre Nations	Hertistrasse 2
5, rue Edouard Belin	CH-8304 Wallisellen
F-44360 Vigneux de Bretagne	

BAXTER is a trademark of Baxter International Inc. TISSEEL is a trademark of Baxter AG. International Inc. TRICOS is a trademark of Baxter Healthcare S.A.

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